

## MAPTrix™-L Laminin mimetic mussel adhesive protein

### Overview

Laminin is a major component of basement membranes, and has numerous biological activities including promotion of cell adhesion, migration, chemotaxis, growth, and differentiation, including neurite outgrowth<sup>1</sup>.

Laminins consist of one  $\alpha$ , one  $\beta$ , and one  $\gamma$  chain, and each chain in the laminin molecule consists of rod-like, globular, and coiled coil regions<sup>1,2</sup>.

The largest chain is  $\alpha$  chain, which contains the long arm on the C-terminal end and a short arm on the N-terminal end. The C-terminal end of the long arm consists of the LG1–5 domains, which are involved in interactions with cellular receptors such as integrins and dystroglycans. The N-terminal end of the short arm is also capable of binding to integrin receptors<sup>3</sup>.


### Product Description

The major laminin's integrin receptors are at least eight integrins ( $\alpha1\beta1$ ,  $\alpha2\beta2$ ,  $\alpha3\beta1$ ,  $\alpha6\beta1$ ,  $\alpha6\beta4$ ,  $\alpha7\beta1$ ,  $\alpha9\beta1$ ,  $\alpha\upsilon\beta3$ ) can bind to laminins<sup>4</sup>.

MAPTrix™-L provides  $\alpha3\beta1$ ,  $\alpha6\beta1$ ,  $\alpha\upsilon\beta3$ , syndecan and dystroglycan binding peptide motifs, derived primarily from N-terminal and globular domain of laminin  $\alpha1$ ,  $\alpha3$  and  $\alpha5$ .

MAPTrix™-L has been demonstrated to have biological activities similar to those of the naturally occurring whole laminin molecule in primary rat hippocampus culture.

**Figure1:** Primary nerve cell cultured on MAPTrix™-L coated surface in serum free conditions



Rat hippocampal cells were primary cultured on MAPTrix™-L and Gibco™ Mouse laminin coated 24-well plates under serum free conditions for 14 days. The MAPTrix™-L showed comparable efficacy for pyramidal dendrite formation of primary nerve cells as compared with mouse-derived natural laminin.

### Characteristics

MAPTrix™-L is produced in Kollodis' proprietary *E.coli* expression system and purified using an ISO compliant manufacturing process.

#### Molecular Weight:

- ~24,000 dalton

#### Formula:

- The product is supplied as a 0.2 mg/mL, 0.5 mg/mL (2.5mg or 5.0mg vial) or 1.0mg/mL aqueous solution in pure water.
- Lyophilized powder is also available upon request

#### Solubility:

- Soluble in a variety of buffers, including water, under a wide range of pH conditions (pH=2~9.0)
- Note: Buffers of media containing  $\text{Ca}^{2+}$  or  $\text{Mg}^{2+}$  added to MAPTrix™ may result in the formation of insoluble aggregates. This will not occur if the buffering capacity of the diluent brings the pH to 9.0 or lower.

### Quality Control

- |                 |  |
|-----------------|--|
| · Purity        | 93% by SDS PAGE  |
| · pH            | 6.0 ~ 7.5  |
| · Endotoxin     | Less than 20 EU/mL per LAL assay.  |
| · Sterility     | Tested and found negative for the presence of bacteria, fungi and mycoplasma                                 |
| · Functionality | The biological activity of laminin peptide is determined in a cell culture assay under serum free conditions |

#### Coating Procedure:

- Transfer desired volume of MAPTrix™-L solution from the vial to a dilution vessel as required.
- Dilute to desired concentration using sodium bicarbonate buffer solution ( $\text{NaHCO}_3$ : 500mM at final concentration) for uniform & even coated surface. A recommended working concentration is 0.1mg/mL. (*Note: Use the recommendation as guidelines to determine the optimal coating conditions for your culture system.*)
- Add appropriate amount of diluted MAPTrix™-L solution to the culture surface
- Incubate at room temperature or 37°C, covered, for 1-3 hours. Best uniform coated surface with 1-2 hr incubation.
- Rinse the coated surfaces carefully with sterile medium or PBS. Avoid scratching the coated surface.
- Refer to the Standard Coating Protocol for details, which can be downloaded at [www.kollodis.com](http://www.kollodis.com)

## Products

Cat. No	Peptide Motif	Receptor	Cat. No	Peptide Motif	Receptor
162041~4	RQVFQVAYIIIIKA (α1 chain)	αvβ3, α5β1	163691~4	GIIFFL (α5 chain)	unknown
162241~4	IKVAV (α1 chain)	αvβ3, α3β1, α6β1,	164111~4	RYVVLPR (β1 chain)	heparin
162261~4	NRWHSIYITRFG (α1 chain)	α6β1	164141~4	YIGSR (β1 chain)	α1β1, α3β1
162291~4	TWYKIAFQRNRK (α1 chain)	α6β1	164211~4	LGTIPG (β1 chain)	67 kDa protein
162321~4	RKRLQVQLSIRT (α1 chain)	syndecan	164421~4	KAFDITYVRLKF (γ1 chain)	αvβ3
162931~4	KNSFMALYLSKGRLVFALG (α3)	syndecan	164601~4	RNIAEIIKDI (γ1 chain)	unknown

For our entire laminin mimetic product listing, please visit our website at [www.kollodis.com](http://www.kollodis.com)

## Storage Conditions:

- Stable for a minimum of 6 months from day of shipment when stored at 2-8°C
- Remaining, unused solution of MAPTrix™ ECM can be stored at 2-8°C with appropriate sealing for 6 months. **DO NOT FREEZE** the remaining solution. However, the remaining material is recommended to be used within 1 month after the vial has been opened.

## References

1. Miner JH, et al., Laminin functions in tissue morphogenesis. *Annu Rev Cell Dev Biol.* 2004; 20:255-84.
2. Colognato H, Form and function: the laminin family of heterotrimers. *Dev Dyn.* 2000; 218(2). p213-34.
3. Madeleine Durbeej, Laminins Cell and Tissue Research (2010) 339(1), 259-268
4. Kim BJ, et al., Mussel adhesive protein fused with cell adhesion recognition motif triggers integrin-mediated adhesion and signaling. *J Biomed Mater Res A.* (2010) 94(3):886-92

## Ordering Information

### USA & Worldwide

AMS Biotechnology

- [www.amsbio.com](http://www.amsbio.com)
- [info@amsbio.com](mailto:info@amsbio.com)
- +1.949.765.8365

Gentaur

- [www.gentaur.com](http://www.gentaur.com)
- [sales@genprice.com](mailto:sales@genprice.com)
- +1.408.472.2934

Kollodis BioSciences

- [www.kollodis.com](http://www.kollodis.com)
- [orders@kollodis.com](mailto:orders@kollodis.com)
- +1.617.283.2182

Sigma-Aldrich

- [www.sigmaaldrich.com](http://www.sigmaaldrich.com)
- +800.325.3010 (within USA)

### Europe

Spain :

- Antibody Bcn
- [www.antibodybcn.com](http://www.antibodybcn.com)
- [info@antibodybcn.com](mailto:info@antibodybcn.com)
- +34.902.220.246

U.K. & The rest :

- AMS Biotechnology
- [www.amsbio.com](http://www.amsbio.com)
- [info@amsbio.com](mailto:info@amsbio.com)
- +44 (0) 1235.232100

Gentaur

- [www.gentaur.com](http://www.gentaur.com)
- [info@gentaur.com](mailto:info@gentaur.com)
- +32.1658.9045

### Asia

China :

- Dakewe Biotech
- [www.dakewe.com](http://www.dakewe.com)
- [info@dakewe.com](mailto:info@dakewe.com)
- +86.755.26650164

4A Biotech Co. Ltd.

- [www.4abio.com](http://www.4abio.com)
- [info@4abio.com](mailto:info@4abio.com)
- +86.400.7060.959

Indonesia : Precision Tech

- [www.pretech.com.sg](http://www.pretech.com.sg)
- [scitech@pretech.com.sg](mailto:scitech@pretech.com.sg)
- +65. 6273.4573

Japan :

Funakoshi Co. Ltd

- [www.funakoshi.co.jp](http://www.funakoshi.co.jp)
- [reagent@funakoshi.co.jp](mailto:reagent@funakoshi.co.jp)
- +81.3.5684.1620

Nacalai Tesque

- [www.nacalai.co.jp](http://www.nacalai.co.jp)
- [info-tech@nacalai.co.jp](mailto:info-tech@nacalai.co.jp)
- +81.75.211.2703

Korea : KDR

- [www.kdr.co.kr](http://www.kdr.co.kr)
- [kdrbio@kdr.co.kr](mailto:kdrbio@kdr.co.kr)
- +82.2.3427.6000

Malaysia/Singapore: Precision Tech

- [www.pretech.com.sg](http://www.pretech.com.sg)
- [scitech@pretech.com.sg](mailto:scitech@pretech.com.sg)
- +65. 6273.4573

Taiwan: Bertec Enterprise Co., Ltd.

- [www.bertec.com.tw](http://www.bertec.com.tw)
- [bertec@bertec.com.tw](mailto:bertec@bertec.com.tw)
- +886.2.2228.1324

For volume ordering or bulk pricing, please contact Kollodis BioSciences or your local distributor.